

Soule Medical™

"Committed to Patient Comfort and Protection"



shown:
SRLE-CC
SUEBS-CC
SRLE-CC
SCB7-CC
(other sizes available)



antimicrobial
LATEX-FREE

ORTHOPEDIC POSITIONERS WITH **microarmor™** **TECHNOLOGY**

ANTIMICROBIAL • FLUID-PROOF • SEAMLESS • DURABLE • RADIOLOGENT • EASY TO CLEAN • REUSABLE



Soule Medical™
Phone: 800.999.2928 • www.soulemedical.com

ORTHOPEDIC POSITIONERS WITH **microarmor™** TECHNOLOGY

Antimicrobial • Fluid-Proof • Seamless • Durable • Radiolucent • Easy to clean • Reusable

• **Microarmor Coating™** is an advanced, durable, vinyl-based coating that encapsulates the foam, forming a high strength bond. Through our advanced application process, this coating provides a robust, fluid-proof barrier without seams or other bacteria harboring areas.

TABLE OF CONTENTS

LOWER EXTREMITY

Leg Ramp Elevator	Page 2
Prone Positioning Wedge	Page 2
Leg Tunnels	Page 2
Acute Positioning Wedge	Page 3
Obtuse Positioning Wedge	Page 3
Knee Bolster	Page 3
Leg Elevator	Page 4
Cylindrical Positioning Rolls	Page 4
Lateral Wedge	Page 4
Podiatry Pad	Page 5
Leg Elevator	Page 5
Well Leg Holder	Page 5

UPPER EXTREMITY

Extremity Wedges	Page 6
Lateral Upper Extremity Positioner	Page 6
Axillary Cradle	Page 6

HEAD RESTS

Head Positioner	Page 7
Donut Head Rest	Page 7

WEDGES

Lateral Positioning Wedge	Page 7
Angular Bolster	Page 8
80° Positioning Wedge	Page 8
40° Positioning Wedge	Page 8

ORTHOPEDIC POSITIONERS WITH **microarmor™** TECHNOLOGY

Leg Ramp Elevator



SRLE-CC	Leg Ramp Elevator	each
SRLE-PP-CC	Ramp/Leg Elevator (Pediatric)	each
SRLE-CC-4	Low Profile Ramp/Leg Elevator	each
SRLE-PP-CC-4	4" Ramp/Leg Elevator (Pediatric)	each

The Leg Ramp Elevator allows you to raise and provide a stable platform for the affected leg.



Prone Positioning Wedge




SRPP-CC	Prone Positioning Wedge	each
----------------	-------------------------	------

Easily elevate the patient's leg while having them in the prone position, providing a working platform for lateral C-Arm imaging.



Leg Tunnels



SLLT-CC	Leg Tunnel (Left)	each
SRLT-CC	Leg Tunnel (Right)	each
SULT-CC	Universal Leg Tunnel	each

Used to provide a stable work platform for operating when the patient is in the lateral position. Contoured tunnel allows the unaffected leg to lay naturally.



ORTHOPEDIC POSITIONERS WITH **microarmor™** TECHNOLOGY

Acute Positioning Wedge



SAWS-CC	Acute Positioning Wedge (Small)	each
SAWM-CC	Acute Positioning Wedge (Medium)	each
SAWL-CC	Acute Positioning Wedge (Large)	each

Holds the knee at acute angle. Gentle curve supports and stabilizes the leg. Eliminates the need for rolled towel bump.



Obtuse Positioning Wedge



SOWS-CC	Obtuse Positioning Wedge (Small)	each
SOWL-CC	Obtuse Positioning Wedge (Large)	each

Holds knee at an obtuse angle. Gentle curve supports and stabilizes the leg. Can assist with fracture reduction. Eliminates the need for rolled towel bump. Small and Large sizes available.



Knee Bolster



SKB-CC	Knee Bolster	each
---------------	--------------	------

Provides gentle elevation and support to the knee.



ORTHOPEDIC POSITIONERS WITH **microarmor™** TECHNOLOGY

Leg Elevator



SERLE-R-CC	Post-op Leg Elevator (Regular)	each
SERLE-R-DF	Post-op Leg Elevator (Regular-Disposable)	box of 3
SERLE-L-CC	Post-op Leg Elevator (Large)	each
SERLE-L-DF	Post-op Leg Elevator (Large-Disposable)	box of 3

Elevates and cradles the leg. Can be used post op to reduce swelling. May also be used to reduce edema.



Cylindrical Positioning Rolls



SCB5-CC	Positioning Bump (5" Diameter)	each
SCB7-CC	Positioning Bump (7" Diameter)	each
SCB9-CC	Positioning Bump (9" Diameter)	each

Provides fast and flexible positioning options for upper or lower extremities. Can be used in both prone or supine positioning. Comes in 5", 7", and 9" diameters.



Lateral Wedge



SLW-R-C	Lateral Wedge (Regular)	each
SLW-S-CC	Lateral Wedge (Small)	each

Provides support and comfort for lateral positioning in the operating room. Ideal for keeping bed-ridden patients repositioned.



ORTHOPEDIC POSITIONERS WITH **microarmor™** TECHNOLOGY

Podiatry Pad



SEH-CC

Podiatry Pad / Extremity Holder

each

Provides stability and elevation of the foot and ankle during podiatry procedures. Can be used to hold upper and lower extremity while prepping surgical site.



Leg Elevator



LP-200-CC

Leg Elevator

each

Elevates and the cradles the leg. Great to use post operative to reduce swelling. May also be used to reduce edema.



Well Leg Holder



WLH-101-CC

Well Leg Holder


each

Supports the non-affected leg during arthroscopic knee procedures.



ORTHOPEDIC POSITIONERS WITH **microarmor™** TECHNOLOGY

Extremity Wedges



SLUEWS-CC	Small Upper Extremity Wedge (Left)	each
SRUEWS-CC	Small Upper Extremity Wedge (Right)	each
SLUEWL-CC	Large Upper Extremity Wedge (Left)	each
SRUEWL-CC	Large Upper Extremity Wedge (Right)	each
SUUELW-CC	Universal Upper Extremity Wedge (Lateral)	each

Provides stable platform for upper extremity when in lateral position. Off-set tunnel allows for C-Arm Imaging and IV access. Angled edge allows space for endotracheal tube. Small and Large sizes available.



Lateral Upper Extremity Positioner



SUEB-S-CC	Upper Extremity - Small	each
SUEB-L-CC	Upper Extremity - Large	each

Allows for easy exposure and provides a stable working area for upper extremity fractures. Cutout allows space for unaffected arm to be extended out from the torso. Angled profile to accommodate intubation devices. Large and small sizes available.



Axillary Cradle



SAC-CC	Axillary Cradle	each
SAC-DF	Axillary Cradle- Disposable	box of 10

Cradled positioner to support the chest during lateral positioning.



ORTHOPEDIC POSITIONERS WITH **microarmor™** TECHNOLOGY

Head Positioner

HR101-CC

Head Positioner 8"x9"x4.5"

each



Donut Head Rest

HR104-CC

Donut Head Rest 9"

each



Lateral Positioning Wedge

WD-106-CC

Wedge: 45°, 24"x7"x7"

each



Supports the patient for lateral and semi-lateral positioning.



ORTHOPEDIC POSITIONERS WITH **microarmor™** TECHNOLOGY

Angular Bolster



WD-114-CC

Angular Bolster

each

Supports the patient. Used to gently elevate and support the knee.



80° Positioning Wedge



WD-122-CC

Wedge: 80°, 21"x10.5"x7.25"

each

80° Positioning Wedge



40° Positioning Wedge



WD-22129-CC

Wedge: 40°, 22"x12"x9"

each

40° Positioning Wedge



ORTHOPEDIC POSITIONERS WITH **microarmor™** TECHNOLOGY

Antimicrobial • Fluid-Proof • Seamless • Durable • Radiolucent • Easy to clean • Reusable

Custom Positioners

To assist you with ALL of your positioning needs

Whether you have a new product idea or just need special sizing or a color option we do not have listed, we can help. Since we are a manufacturer, we are able to provide custom products without charging a custom price.

Minimum order quantities may apply.





Soule MedicalTM

"Committed to Patient Comfort and Protection"



OR Patient Positioning Systems

microarmorTM

TECHNOLOGY

ANTIMICROBIAL • FLUID-PROOF • SEAMLESS • DURABLE • RADIOLOGENT • EASY TO CLEAN • REUSABLE



Soule MedicalTM
Phone: 800.999.2928 • www.soulemedical.com

